

PACE 101 H

Externer Demand-Herzschrittmacher
External Demand Pacemaker



Externer Demand-Herzschrittmacher	External Demand-Pacemaker	
<p>PACE 101 H gewährleistet eine zuverlässige und effektive temporäre Stimulation bei</p> <ul style="list-style-type: none"> • Reizbildungs- u. Reizleitungsstörungen • Bradykardien nach herzchirurgischem Eingriff • Implantation bzw. Austausch von Herzschrittmachern <p>Besondere Eigenschaften:</p> <ul style="list-style-type: none"> • Einstellung der Parameter per Knopf • Berührungsgeschützte Spannzangen für Stecker mit 0,9 - 2,0 mm Durchmesser • Langlebige Batterie (1000 h) • Batterieüberwachung mit akustischer und optischer Warnanzeige • Grundfrequenz bis zu 180 ppm für pädiatrische Anwendung • Overdrive-Stimulation durch Verdoppelung bzw. Vervierfachung der Stimulationsfrequenz zur Terminierung von atrialen Tachykardien <p>Technische Daten</p>	<p>PACE 101 H provides reliable and effective temporary stimulation for</p> <ul style="list-style-type: none"> • Arrhythmia and cardiac conduction disorders • Bradycardia after cardiac surgery • Pacemaker implant or replacement procedures <p>Special Features:</p> <ul style="list-style-type: none"> • Dial setting of parameters • Protected collets for pins of 0.9 - 2.0 mm diameter • Long-life battery (1000 h) • Battery monitoring with acoustic and optical warning • Basic rate up to 180 ppm for pediatric application • High rate stimulation by doubling or quatrupling the stimulation rate to terminate atrial tachycardia <p>Technical Data</p>	
Modi	Modes	A00, AAI, V00, VVI
Grundfrequenz	Basic rate	30 ... 180 ppm
Overdrive-Stimulation	High rate	async., freq. x2, x4
Stimulationsamplituden	Stimulation amplitude	0.3 ... 12 V
Impulsdauer	Pulse duration	0.75 ms
Empfindlichkeit	Sensitivity	1.0 ... 20 mV, [∞]
Refraktärzeit	Refractory period	250 ms
Batterie Lebensdauer	Battery life span	1000 h
Gewicht mit Batterie	Weight with battery	170 g

Manufacturer: **Osyпка Medical GmbH**
D-12277 Berlin, Telefon +49 (30) 741 60 58

Dr. Osypka GmbH, Medizintechnik
Earl-H.-Wood-Str. 1, 79618 Rheinfelden, Germany
Tel. +49 (7623) 74 05-0, Fax +49 (7623) 74 05-213

CE 0123

www.osypka.de